

PYROLYSIS

CAPILLARY GAS CHROMATOGRAPHY ANALYSIS OF SOME MODERN PAINT FINISHES

INTRODUCTION

Capillary gas chromatographic profiles of three commercially available paint finishes were obtained by a simple and rapid pyrolysis technique.


The samples run were a water based flat white plastic paint (household use), a white, gloss enamel (all purpose style) and a blue, acrylic paint (used in automotive industry). Sample preparation was simple using an SGE sample "Pelletiser P-1" and a "Solids Injector" system.

Results exhibit excellent reproducibility and extremely high sensitivity.

All samples were run under the following conditions.

INSTRUMENTATION AND CONDITIONS:

Pyrolysis Temp.:	800°C
Column:	25QC3/BP1/0.5, 25m, 0.32mm ID
Stationary Phase:	Crosslinked methyl silicone, BP1
Film thickness:	0.5 micron
Initial Temp.	50°C for 2min
Program Rate:	5°C/min
Final Temp.:	200°C for 10min
Detector:	Flame Ionization
Carrier Gas:	Hydrogen
Linear Carrier Velocity:	35cm/sec @ 50°C
Split Ratio:	100:1
Injection Mode:	Pyrolysis Split



SGE International Pty. Ltd.
Toll Free: 1800 800 167
Tel: +61 (0) 3 9837 4200
Fax: +61 (0) 3 9874 5672
Email: support@sge.com

SGE Europe Ltd
Tel: +44-1908 568 844
Fax: +44-1908 566 790
Email: uk@sge.com

SGE Europe Ltd
Tel: +33 1 6929 8090
Fax: +33 1 6929 0925
Email: france@sge.com

SGE, Incorporated (USA)
Toll Free: (800) 945 6154
Tel: +1-512-837 7190
Fax: +1-512-836 9159
Email: usa@sge.com

SGE GmbH
Tel: +49 6151 860486
Fax: +49 6151 860489
Email: germany@sge.com

SGE Japan Inc.
Tel: +81 45 222 2885
Fax: +81 45 222 2887
Email: japan@sge.com

SGE India
Tel: +91 (22) 2471 5896
Fax: +91 (22) 2471 6592
Email: sgeindia@vsnl.com

SGE China
Tel: +86 (10) 6588 8666
Fax: +86 (10) 6588 6577
Email: sales@jindustries.com.cn

SGE Gulf
Tel: +971-6-557 3341
Fax: +971-6-557 3541
Email: gulfsupport@sge.com