

















MARKET FOCUSED - Environmental PYRETHROID ANALYSIS

Pyrethroids are potent insecticides widely used in agriculture, disease control and household products. Nearly all domestic fly aerosols use pyrethroids as the active ingredient whether these are surface or airborne sprays. Pyrethroids act by interfering with the insect's nervous system and, in high concentrations, affect humans in a similar way.

Synthetic pyrethroids are also known as endocrine disruptors. Other toxicological properties include: liver damage when subjected to long term exposure, allergic reactions and asthmatic attacks. The USEPA has found that many of the pyrethroids may also be carcinogenic. For this reason, the levels of synthetic pyrethroids in foods and the environment are always under close scrutiny. Analysis by GC can be complex because many of the pyrethroids can have up to 8 optical isomers thus more than one peak can represent the same pyrethroid. Therefore, overlapping of signals can be a problem (e.g. Cypermethrin with 3 chiral centers can have up to 8 optical isomers, which include 4 diastereoisomers) (Figure 3). These identification problems have been made a thing of the past with SGE's 5% phenyl BPX5, 35% phenyl BPX35, and 50% phenyl BPX50 columns. Each column gives superb separation of the various isomers of the 16 components (Table 1) in less than 33 minutes, allowing a fast turnaround time of samples. With an upper temperature limit of 360°C on all three columns, high boiling point contaminants can be "baked out" without damaging the stationary phase. This ensures that contaminants won't interfere with subsequent analyses. The bleed profiles at 300°C are excellent and the exceptional peak shape of each pyrethroid indicates a high degree of inertness on all three columns.

Table 1 shows the structures of the various Pyrethroids analyzed

| PYRETHROID STRUCTURES | | | |
|-----------------------|--|-------------------|---|
| Pyrethroid | Compound Structure | Pyrethroid | Compound Structure |
| 1. Natural Pyrethrums |  | 9. Permethrin |  |
| 2. Allenthin |  | 10. Cyfluthrin |  |
| 3. Bioresmethrin |  | 11. Cypermethrin |  |
| 4. Resmethrin |  | 12. Flucythrinate |  |
| 5. Bifenthrin |  | 13. Fenvalerate |  |
| 6. Fenproprathrin |  | 14. Esfenvalerate |  |
| 7. Tetramethrin |  | 15. Tralomethrin |  |
| 8. Sumithrin |  | 16. Deltamethrin |  |

BPX5

5% phenyl (equiv) polysilphenylene-siloxane

BPX5

Replaces

DB-5
DB-5MS
DB-5.625
XTI-5
Rtx-5ms
Ultra-2
HP-5
HP-5MS
HP5-TA
SPB-5
MDN-5S
CP-Sil 8CB
Rtx-5Sil MS
AT-5
CP-Sil 8CB M


Figure 1 shows the separation of 16 pyrethroids on BPX5. Note the superb bleed profile at 300°C. BPX5 gives excellent separation of the synthetic pyrethroids. Breakdown in the injection port of the thermally labile pyrethroids to their parent compounds such as tralomethrin (16) to deltamethrin (17) and esfenvalerate (14) to fenvalerate (15) are the only coelutions seen.

Figure 1

Phase: BPX5, 0.25µm film
Sample: 100ppm in Dichloromethane
Column: 30m x 0.25 mm ID
Initial Temp: 110°C, 1 min.
Rate 1: 25°C/min to 150°C,
Rate 2: 12°C/min to 260°C
Rate 3: 15°C/min to 310°C
Final Temp: 310°C, 6 min.
Detector Type: Mass Spectrometer
Carrier Gas: He, 9.5 psi
Carrier Gas Flow: 0.9 mL/min.
Constant Flow: On
Average Linear Velocity: 35 cm/sec at 110°C
Injection Mode: Split
Split Ratio: 50:1
Injection Volume: 1 mL
Injection Temperature: 250°C
Autosampler: No
Liner Type: 4 mm ID Single Taper Liner
Liner Part Number: 092017
Column Part Number: 054101

COMPONENTS

1. Natural Pyrethrums
2. Allethrin
3. Bioresmethrin
4. Resmethrin
5. Bifenthrin
6. Fenpropathrin
7. Tetramethrin
8. Sumithrin
9. Permethrin
10. Cyfluthrin
11. Cypermethrin
12. Flucythrinate
13. Fenvalerate
14. Esfenvalerate
15. Tralomethrin
16. Deltamethrin



BPX35

35% phenyl (equiv) polysilphenylene-siloxane

BPX35
Replaces



DB-35
DB-35MS
Rtx-35
Rtx-35ms
HP-35
HP-35MS
SPB-35
MDN-35
AT-35

Figure 2 shows excellent separation of the 16 synthetic pyrethroids on BPX35. The bleed at 300°C is excellent and the signal to noise in full scan mode is 60:1 for 10ng on column.


Figure 3 shows a chromatogram of the synthetic pyrethroid cypermethrin. Note the 4 peaks seen here representing the 4 diastereoisomers of cypermethrin. The structure of cypermethrin shown here shows the 3 chiral centers (*) of cypermethrin.

Figure 2 & Figure 3

| | |
|---------------------------------|---------------------------|
| Phase: | BPX35 0.25µm film |
| 16 Pyrethroids: | 10ppm in methanol |
| Column: | 30m x 0.25mm ID |
| Initial Temp: | 50°C, 1 min |
| Rate 1: | 30°C/min to 200°C |
| Rate 2: | 4°C/min to 300°C |
| Final Temp: | 300°C for 5 min |
| Detector Type: | Mass Spectrometer |
| Carrier Gas: | He, 6.5 psi |
| Carrier Gas Flow: | 0.9mL/min |
| Constant Flow: | On |
| Average Linear Velocity: | 35cm/sec at 50°C |
| Injection Mode: | Spiltless |
| Purge on Time: | 0.5min |
| Purge on (Split) Vent: | 60mL/min |
| Injection Volume: | 1µL |
| Injection temperature: | 250°C |
| Autosample: | No |
| Liner Type: | 4mm ID Double Taper Liner |
| Liner Part Number: | 092018 |
| Column part Number: | 054701 |



Refer to Table 1 for Peak Identification



BPX50

50% phenyl (equiv) polysilphenylene-siloxane

Figure 4 show the separation of the synthetic pyrethroids as run on BPX50. Note the run time on this chromatogram is 16 minutes for complete elution of the 16 pyrethroids. This allows analytical laboratories to have higher sample throughput.

BPX50 Replaces


OV-17
SP-2250
DB-17
DB-17MS
DB-17ht
Rtx-50
SPB-50
HP-50+
HP-17
AT-50
007-17

Figure 4

Phase: BPX50 0.25 μ m film
16 Pyrethroids: 10ppm in methanol
Column: 30m x 0.25mm ID
Initial Temp: 50°C, 1 min
Rate 1: 30°C/min to 200°C
Rate 2: 4°C/min to 300°C
Final Temp: 300°C for 5 min
Detector Type: Mass Spectrometer
Carrier Gas: He, 6.5psi
Carrier Gas Flow: 0.9mL/min
Constant Flow: On
Average Linear Velocity: 35cm/sec at 50°C
Injection Mode: Splitless
Purge On Time: 0.5min
Purge on (Split) Vent: 60mL/min
Injection Volume: 1 μ L
Injection Temperature: 250°C
Autosampler: No
Liner Type: 4mm ID Double Taper Liner
Liner Part Number: 092018
Column Part Number: 054751

COMPONENTS

1. Natural Pyrethrums
2. Allethrin
3. Bioresmethrin
4. Resmethrin
5. Bifenthrin
6. Fenpropathrin
7. Tetramethrin
8. Sumithrin
9. Permethrin
10. Cyfluthrin
11. Cypermethrin
12. Flucythrinate
13. Fenvalerate
14. Esfenvalerate
15. Tralomethrin
16. Deltamethrin



ORDERING INFORMATION

| | ID(mm) (mm) | Length (m) | Film Thickness (µm) | Temperature Limits (°C) | Part No. | | |
|------|----------------|---------------|------------------------|----------------------------|----------------|--------------|--------|
| BPX5 | 0.10 | 10 | 0.10 | -40 to 360/370 | 054099 | | |
| | | 12 | 0.25 | -40 to 360/370 | 054112 | | |
| | | 25 | 0.25 | -40 to 360/370 | 054113 | | |
| | 0.22 | 30 | 0.25 | -40 to 360/370 | 054142 | | |
| | | 50 | 0.25 | -40 to 360/370 | 054114 | | |
| | | 15 | 0.25 | -40 to 360/370 | 054100 | | |
| | | 15 | 1.00 | -40 to 360/370 | 054121 | | |
| | | 30 | 0.25 | -40 to 360/370 | 054101 | | |
| | | 30 | 0.50 | -40 to 360/370 | 0541025 | | |
| | 0.25 | 30 | 1.00 | -40 to 360/370 | 054122 | | |
| | | 60 | 0.25 | -40 to 360/370 | 054102 | | |
| | | 15 | 0.25 | -40 to 360/370 | 054144 | | |
| | | 15 | 1.00 | -40 to 360/370 | 054152 | | |
| | | 25 | 0.25 | -40 to 360/370 | 054119 | | |
| | | 25 | 0.50 | -40 to 360/370 | 054125 | | |
| | | 30 | 0.25 | -40 to 360/370 | 054145 | | |
| | | 30 | 0.50 | -40 to 360/370 | 0541205 | | |
| | | 30 | 1.00 | -40 to 360/370 | 054145 | | |
| | | 60 | 0.25 | -40 to 360/370 | 054146 | | |
| | | 60 | 1.00 | -40 to 360/370 | 054154 | | |
| | | 0.53 | 12 | 1.00 | -40 to 360/370 | 054130 | |
| | 15 | | 1.00 | -40 to 360/370 | 054147 | | |
| | 15 | | 1.50 | -40 to 350/360 | 0541347 | | |
| | 15 | | 3.00 | -40 to 350/360 | 054159 | | |
| | 25 | | 1.00 | -40 to 360/370 | 054131 | | |
| | 30 | | 0.50 | -40 to 360/370 | 0541345 | | |
| | 30 | | 1.00 | -40 to 360/370 | 054148 | | |
| | 30 | | 1.50 | -40 to 350/360 | 0541348 | | |
| | 30 | | 3.00 | -40 to 350/360 | 054160 | | |
| | 60 | | 1.00 | -40 to 360/370 | 054158 | | |
| | 60 | | 3.00 | -40 to 350/360 | 054164 | | |
| | BPX35 | | 0.10 | 10 | 0.1 | 0 to 360/370 | 054699 |
| | | | | 15 | 0.25 | 0 to 360/370 | 054713 |
| | | | | 30 | 0.25 | 0 to 360/370 | 054714 |
| | | | 0.22 | 15 | 0.25 | 0 to 360/370 | 054700 |
| | | | | 15 | 1 | 0 to 360/370 | 054703 |
| 30 | | | | 0.25 | 0 to 360/370 | 054701 | |
| 30 | | | | 0.5 | 0 to 360/370 | 0547025 | |
| 30 | | 1 | | 0 to 360/370 | 054704 | | |
| 60 | | 0.25 | | 0 to 360/370 | 054702 | | |
| 0.25 | | 60 | 1 | 0 to 360/370 | 054705 | | |
| | | 15 | 0.25 | 0 to 360/370 | 054723 | | |
| | | 15 | 0.5 | 0 to 360/370 | 054718 | | |
| | | 15 | 1 | 0 to 360/370 | 054716 | | |
| | | 30 | 0.25 | 0 to 360/370 | 054724 | | |
| | | 30 | 0.5 | 0 to 360/370 | 0547158 | | |
| | | 30 | 1 | 0 to 360/370 | 054717 | | |
| | | 60 | 0.25 | 0 to 360/370 | 054725 | | |
| | | 0.32 | 15 | 0.5 | 0 to 360/370 | 054734 | |
| | 15 | | 1 | 0 to 360/370 | 054736 | | |
| | 30 | | 0.5 | 0 to 360/370 | 054735 | | |
| | 30 | | 1 | 0 to 360/370 | 054737 | | |
| 0.53 | 10 | | 0.10 | 0 to 360/370 | 054740 | | |
| | 15 | | 0.25 | 0 to 360/370 | 054750 | | |
| | 30 | 0.25 | 0 to 360/370 | 054751 | | | |
| | 60 | 0.25 | 0 to 360/370 | 054752 | | | |
| | 0.32 | 15 | 0.25 | 0 to 360/370 | 054760 | | |
| | | 30 | 0.25 | 0 to 360/370 | 054761 | | |
| 0.53 | 15 | 0.50 | 0 to 360/370 | 054770 | | | |
| | 30 | 0.50 | 0 to 360/370 | 054771 | | | |
| | 30 | 1.00 | 0 to 360/370 | 054772 | | | |
| | SolGel-1ms™ | 30 | 0.25 | 0 to 370/380 | 054795 | | |
| | | 60 | 0.25 | 0 to 370/380 | 054793 | | |
| | | 30 | 0.25 | 0 to 370/380 | 054798 | | |
| | 0.32 | 60 | 0.25 | 0 to 370/380 | 054794 | | |

SUMMARY:

BPX5, BPX35 and BPX50 columns show unparalleled performance for the separation of synthetic pyrethroids. These columns can be conditioned at the end of each analysis to remove any high boiling point contaminants without any degradation to the stationary phase and give excellent signal-to-noise with very little on-column breakdown. BPX5, BPX35 and BPX50 are the columns of choice for all pyrethroid analyses.

For information on SGE's complete range of GC accessories, please refer to our latest catalog, or contact your nearest SGE office or distributor.



Publication No. TA-0076-C Rev: 02 05/01

| | | | | |
|--|--|---|---|--|
| SGE International Pty. Ltd. Tel: +61 (0) 3 9837 4200 Fax: +61 (0) 3 9874 5672 email: support@sge.com | SGE Incorporated (USA) Toll Free: 800 945 6154 Fax: +1 512 836 9159 email: usa@sge.com | SGE France Sarl Tel: +33 (0) 1 6929 8090 Fax: +33 (0) 1 6929 0925 email: france@sge.com | SGE Deutschland GmbH Tel: +49 (0) 6151 860486 Fax: +49 (0) 6151 860489 email: germany@sge.com | SGE Japan Inc. Tel: +81 (45) 222 2885 Fax: +81 (45) 222 2887 email: japan@sge.com |
| SGE Australia Pty. Ltd. Toll Free: 1800 800 167 Fax: (03) 9874 5672 email: sales@sge.com | SGE Europe Ltd. (UK) Tel: +44 (0) 1908 568 844 Fax: +44 (0) 1908 566 790 email: uk@sge.com | SGE Italia Srl. Tel: +39 06 4429 0206 Fax: +39 06 4429 0724 email: sge.italia@tin.it | SGE China Service Centre Tel: +86 (10) 6588 8666 Fax: +86 (10) 6588 6577 email: sge@bj.col.com.cn | SGE India Tel: +91 (022) 471 5896 Fax: +91 (022) 471 6592 email: sgeindia@bom4.vsnl.net.in |