

UNIVERSAL SYRINGES™

...MORE THAN JUST FLEXIBLE

The Universal syringe features a Super-elastic alloy plunger and a cone tip needle which extends trouble-free syringe life.

Universal syringes have gained wide acceptance and are the perfect upgrade for standard syringe types. The indestructible super-elastic plunger returns to its original shape if twisted or bent. The cone tip needles offer extended septum life through reduced coring during multiple injections.

Universal syringes are available in 5µL and 10µL volumes. The 10µL syringe comes in two needle versions, 0.63mm OD and 0.47mm OD. The 0.63mm oD needle is recommended for applications where there is high probability of needle bending, such as student use.



PRODUCT DATA

ORDERING INFORMATION

Volume	Length (mm)	Gauge	OD (mm)	Tip Style	Syringe Code	Part No.
Fixed Needle						
5µL	50	26	0.47	Bevel	5FX	001100
5µL	50	26	0.47	Cone	5FX-5C	001105
10µL	50	26	0.47	Bevel	10FX	002100
10µL (6 pack)	50	26	0.47	Bevel		002130
10µL	50	26	0.47	Cone	10FX-5/0.63C	002108
10µL	70	26	0.47	Bevel	10FX-7	002101
Removable Needle						
5µL	50	26	0.47	Bevel	5RX	001150
Spare Needles (packet of 5)						
10µL	50	26	0.47	Bevel	10RX	002150
10µL (6 pack)	50	26	0.47	Bevel		002180
Spare Needles (packet of 5)						037110

Universal Syringe™: Trademark of SGE International

Publication No. PD-0062-S Rev: 02 03/02



SGE International Pty. Ltd. Toll Free: 1800 800 167 Tel: +61 (0) 3 9837 4200 Fax: +61 (0) 3 9874 5672 email: support@sge.com	SGE Incorporated (USA) Toll Free: 800 945 6154 Fax: +1 512 836 9159 email: usa@sge.com	SGE (France) Sarl Tel: +33 (0) 1 6929 8090 Fax: +33 (0) 1 6929 0925 email: france@sge.com	SGE Deutschland GmbH Tel: +49 (0) 6151 860486 Fax: +49 (0) 6151 860489 email: germany@sge.com	SGE Japan Inc. Tel: +81 (45) 222 2885 Fax: +81 (45) 222 2887 email: japan@sge.com
	SGE Europe Ltd. (UK) Tel: +44 (0) 1908 568 844 Fax: +44 (0) 1908 566 790 email: uk@sge.com	SGE (Italia) Srl. Tel: +39 06 4429 0206 Fax: +39 06 4429 0724 email: sge.italia@tin.it	SGE China Service Centre Tel: +86 (10) 6588 8666 Fax: +86 (10) 6588 6577	SGE India Tel: +91 (022) 471 5896 Fax: +91 (022) 471 6592 email: sgeindia@vsnl.net.in