


Syringe Multipacks

Safety in Numbers

When you purchase SGE's Syringe Multipacks, you **SAVE MONEY** and reduce your cost per analysis without compromising on quality.



6 pack Syringes



10 pack Syringes

UNIQUE ADVANTAGES

- **HIGHLY ECONOMICAL**
SGE offers a **10% saving** on 6 packs and a **25% saving** on 10 packs.
- **CONVENIENT**
Replacements are always at your fingertips so there is never a delay in results.
- **ACCURACY**
Syringes are accurate to +/-1% of dispensed volume.
- **OPTIONS**
 - Plunger Protection
 - Autosampler
 - SuperfleX - Ti-Ni plunger
 - Gas Tight
- **Available for manual or autosampler injection.**
- **LONG LIFE SYRINGES**
SGE's unique plunger options prevent plunger bending and seizing.

APPLICATIONS

SGE Syringe Multipacks are ideal for all GC (manual or autosampler) and HPLC laboratories, especially those involved in high turnover chromatography.

SYRINGE SPECIFICATIONS

Maximum Pressure	1,000psi	
Temperature	Fixed Needle Syringe	5-70°C
	Removable Needle Syringe	5-120°C
Accuracy of dispensed Volume	+/- 1%	
Sterilisation Mode	Autoclave (remove Teflon tip plungers from barrel) Gas (ethylene oxide) Disinfectants	

Syringe Ordering Information

SYRINGES FOR MANUAL INJECTION - WITH PLUNGER PROTECTION

Part No.	Pack Size	Description	Needle OD	Needle Length	Tip Style
002030	Six Pack	10µl fixed needle	0.47mm	50mm	Bevel
002033	Ten Pack	10µl fixed needle	0.47mm	50mm	Bevel
002080	Six Pack	10µL removable needle	0.47mm	50mm	Bevel
002202	Six Pack	10µL Gas Tight fixed needle	0.47mm	50mm	Bevel
002252	Six Pack	10µL Gas Tight removable needle	0.47mm	50mm	Bevel

SYRINGES FOR AUTOSAMPLERS

(Suitable for Agilent Technologies (HP) & Carlo Erba/Fisons Autosamplers)

Part No.	Pack Size	Description	Needle OD	Needle Length	Tip Style
001804	Six Pack	5µL HP7673 fixed needle	0.47mm	42mm	Cone
001814	Six Pack	5µL HP7673 fixed needle	0.63mm	42mm	Cone
002804	Six Pack	10µL HP7673 fixed needle	0.47mm	42mm	Cone
002814	Six Pack	10µL HP7673 fixed needle	0.63mm	42mm	Cone
002986	Six Pack	10µL CTC/Fisons fixed needle	0.47mm	50mm	Cone

SYRINGES FOR HPLC INJECTION - WITH PLUNGER PROTECTION

(Suitable for Rheodyne & Valco Valves)

Part No.	Pack Size	Description	Needle OD	Needle Length	Tip Style
002315	Six Pack	10µL LC fixed needle	0.028"	2"	LC

Mentioned trademarks are registered trademarks of their respective owners. Specifications are subject to change without notice.



© Copyright 2001. SGE International Pty. Ltd. All rights reserved. Publication No. PD-0117-S Rev:01 01/01

SGE International Pty. Ltd. Toll Free: 1800 800 167 Tel: +61 (0) 3 9837 4200 Fax: +61 (0) 3 9874 5672 email: support@sge.com	SGE Incorporated (USA) Toll Free: 800 945 6154 Fax: +1 512 836 9159 email: usa@sge.com	SGE (France) Sarl Tel: +33 (0) 1 6929 8090 Fax: +33 (0) 1 6929 0925 email: france@sge.com	SGE Deutschland GmbH Tel: +49 (0) 6151 860486 Fax: +49 (0) 6151 860489 email: germany@sge.com	SGE Japan Inc. Tel: +81 (45) 222 2885 Fax: +81 (45) 222 2887 email: japan@sge.com.au
	SGE Europe Ltd. (UK) Tel: +44 (0) 1908 568 844 Fax: +44 (0) 1908 566 790 email: uk@sge.com	SGE (Italia) Srl. Tel: +39 06 4429 0206 Fax: +39 06 4429 0724 email: sge.italia@tin.it	SGE China Service Centre Tel: +86 (10) 6588 8666 Fax: +86 (10) 6588 6577 email: sge@bj.col.com.cn	SGE India Tel: +91 (022) 471 5896 Fax: +91 (022) 471 6592 email: sgeindia@vsnl.com