

Injection Port | Inlet Liner Double Taper

Spot the Difference in the SGE Lineup!

Part Number	Part Description and Detail
Agilent	
09201801	Agilent Inlet Liner 4 mm ID Double Taper, Length 78.5 mm PK1
092018	Agilent Inlet Liners 4 mm ID Double Taper, Length 78.5 mm PK5
092230	Agilent Inlet Liners 4 mm ID Double Taper, Length 78.5 mm PK25

Note - Bruker/Varian 1077 injectors and Thermo Trace 1300 series use 78.5 mm "Agilent" Liners

For more information about this product visit www.sge.com or contact techsupport@sge.com.