

# Syringes | Guided Plunger

## Diamond Syringe - Syringe Brilliance


Guided plunger syringes are the most rugged syringe available, making them perfect for industrial environments.

### Features and benefits

- Extended barrel guides plunger during injection.
- Robust and rugged.

### Recommended applications

Ideal for student and industrial use.



### Product specifications

Available in 5  $\mu$ L to 10  $\mu$ L capacities.

Accuracy and Reproducibility	$\pm 1$ % (dispensed volume)
Temperature Range	5-70 $^{\circ}$ C Fixed Needle, 5-120 $^{\circ}$ C Removable Needle
International Standards Traceability	

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)