

Injection Port | Inlet Liner ConnecTite

Spot the Difference in the SGE Lineup!

Part Number	Part Description and Detail
Agilent	
092324	ConnecTite Agilent Inlet Liners Length 78.5 mm PK5
092325	ConnecTite Agilent Inlet Liners with Top Hole, Length 78.5 mm PK5
092326	ConnecTite Agilent Inlet Liners with Bottom Hole, Length 78.5 mm PK5
Bruker/Varian	
09233901	ConnecTite Bruker/Varian 1075/1077 Inlet Liner Length 72 mm PK1
092339	ConnecTite Bruker/Varian 1075/1077 Inlet Liners Length 72 mm PK5
09234001	ConnecTite Bruker/Varian 1075/1077 Inlet Liner with Top Hole, Length 72 mm PK1
092340	ConnecTite Bruker/Varian 1075/1077 Inlet Liners with Top Hole, Length 72 mm PK5
09234101	ConnecTite Bruker/Varian 1075/1077 Inlet Liner with Bottom Hole, Length 72 mm PK1
092341	ConnecTite Bruker/Varian 1075/1077 Inlet Liners with Bottom Hole, Length 72 mm PK5
09203001	ConnecTite Bruker/Varian 1093/1094 Inlet Liner 0.8 mm ID, Length 54 mm (Rest 0.25 mm) PK1
092030	ConnecTite Bruker/Varian 1093/1094 Inlet Liners 0.8 mm ID, Length 54 mm (Rest 0.25 mm) PK5
09203401	ConnecTite Bruker/Varian 1093/1094 Inlet Liner 0.8 mm ID, Length 54 mm (Rest 0.5 mm) PK1
092034	ConnecTite Bruker/Varian 1093/1094 Inlet Liners 0.8 mm ID, Length 54 mm (Rest 0.5 mm) PK5
PerkinElmer	
092344	ConnecTite PerkinElmer Inlet Liners Length 92 mm PK5
092345	ConnecTite PerkinElmer Inlet Liners with Top Hole, Length 92 mm PK5
09234601	ConnecTite PerkinElmer Inlet Liner with Bottom Hole, Length 92 mm PK1
092346	ConnecTite PerkinElmer Inlet Liners with Bottom Hole, Length 92 mm PK5
Shimadzu	
09232901	ConnecTite Shimadzu Inlet Liner Length 95 mm PK1
092329	ConnecTite Shimadzu Inlet Liners Length 95 mm PK5
09233001	ConnecTite Shimadzu Inlet Liner with Top Hole, Length 95 mm PK1
092330	ConnecTite Shimadzu Inlet Liners with Top Hole, Length 95 mm PK5
09233101	ConnecTite Shimadzu Inlet Liner with Bottom Hole, Length 95 mm PK1
092331	ConnecTite Shimadzu Inlet Liners with Bottom Hole, Length 95 mm PK5
Thermo Scientific	
09233401	ConnecTite Thermo Scientific Inlet Liner Length 105 mm PK1
092334	ConnecTite Thermo Scientific Inlet Liners Length 105 mm PK5
09233501	ConnecTite Thermo Scientific Inlet Liner with Top Hole, Length 105 mm PK1
092335	ConnecTite Thermo Scientific Inlet Liners with Top Hole, Length 105 mm PK5
09233601	ConnecTite Thermo Scientific Inlet Liner with Bottom Hole, Length 105 mm PK1
092336	ConnecTite Thermo Scientific Inlet Liners with Bottom Hole, Length 105 mm PK5

Note - Bruker/Varian 1077 injectors and Thermo Trace 1300 series are using 78.5 mm "Agilent" Liners

For more information about this product visit www.sge.com or contact techsupport@sge.com.