

Connections | HPLC Unions

Perfect Connections, Every Time

Part Number	Part Description and Detail
ProteCol Stainless Steel Unions	
211010	1/16 " to 1/16 " (1 Union Only - No Ferrules) ProteCol Stainless Steel Unions
212006	1/16 " to 1/16 " (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
212007	1/16 " to 1/32 " (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
212008	1/16 " to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
222005	1/32 " to 1/32 " (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
222006	1/32 " to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
222007	0.36 mm to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol Stainless Steel Unions
ProteCol PEEK Unions	
212009	1/16 " to 1/16 " (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
212010	1/16 " to 1/32 " (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
212011	1/16 " to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
222009	1/32 " to 1/32 " (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
222010	1/32 " to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
222011	0.36 mm to 0.36 mm (2 Unions & 4 PEEK Ferrules) ProteCol PEEK Unions
ProteCol PEEK Ferrules	
213006	1/16 " to 1/16 ", ProteCol PEEK Ferrules PK5
213007	1/16 " to 1/32 ", ProteCol PEEK Ferrules PK5
213008	1/16 " to 0.36 mm, ProteCol PEEK Ferrules PK5
223005	1/32 " to 1/32 ", ProteCol PEEK Ferrules PK5
223006	1/32 " to 0.36 mm, ProteCol PEEK Ferrules PK5
223007	0.36 mm to 0.36 mm, ProteCol PEEK Ferrules PK5

For more information about this product visit www.sge.com or contact techsupport@sge.com.