

## Connections | GC Ferrules Graphite Vespel®

Part Number	Part Description and Detail
072671	15% Graphite/85% Vespel Ferrule with 0.8 mm Hole for 1/8 " Fitting, PK10
072673	15% Graphite/85% Vespel Ferrule with 1/8 " Hole for 1/16 " Fitting, PK10
072674	15% Graphite/85% Vespel Ferrule with 2 x 0.8 mm Holes for 1/8 " Fitting, PK10
072694	15% Graphite/85% Vespel Mini Union Ferrule with No Hole, PK10
072695	15% Graphite/85% Vespel Mini Union Ferrule with 0.3 mm Hole, PK10
072696	15% Graphite/85% Vespel Mini Union Ferrule with 0.4 mm Hole, PK10
072697	15% Graphite/85% Vespel Mini Union Ferrule with 0.5 mm Hole, PK10
072698	15% Graphite/85% Vespel Mini Union Ferrule with 0.8 mm Hole, PK10
073105	15% Graphite/85% Vespel Short Ferrule with No Hole, for Agilent Injectors, PK10
073107	15% Graphite/85% Vespel Short Ferrule with 0.3 mm Hole, for Agilent Injectors, PK10
073109	15% Graphite/85% Vespel Short Ferrule with 0.4 mm Hole, for Agilent Injectors, PK10
073111	15% Graphite/85% Vespel Short Ferrule with 0.5 mm Hole, for Agilent Injectors, PK10
073113	15% Graphite/85% Vespel Short Ferrule with 0.8 mm Hole, for Agilent Injectors, PK10
0726547	15% Graphite/85% Vespel Ferrule with 0.8 mm Hole, for Shimadzu QP5000-II & QP5050 Standard MS Interface, PK10
0726548	15% Graphite/85% Vespel Ferrule with 0.7 mm Hole for Carlo Erba Instruments, Not Suitable for 8000 Series and Trace/Focus GC, PK10
0726549	15% Graphite/85% Vespel Ferrule with 0.4 mm Hole for Carlo Erba Instruments, Not Suitable for 8000 Series and Trace/Focus GC, PK10
0726557	15% Graphite/85% Vespel Ferrule with 0.5 mm Hole for Carlo Erba Instruments, Not Suitable for 8000 Series and Trace/Focus GC, PK10
0726560	15% Graphite/85% Vespel Ferrule with 0.3 mm Hole, for Shimadzu QP5000-II & QP5050 Standard MS Interface, PK10
0726561	15% Graphite/85% Vespel Ferrule with 0.4 mm Hole, for Shimadzu QP5000-II & QP5050 Standard MS Interface, PK10
0726562	15% Graphite/85% Vespel Ferrule with 0.5 mm Hole, for Shimadzu QP5000-II & QP5050 Standard MS Interface, PK10
0726563	15% Graphite/85% Vespel Ferrule with 0.4 mm Hole, for Shimadzu QP5000-I Standard MS Interface, PK10
0726564	15% Graphite/85% Vespel Ferrule with 0.5 mm Hole, for Shimadzu QP5000-I Standard MS Interface, PK10
0726591	15% Graphite/85% Vespel Vg/Micromass Ferrule with 0.3/0.4 mm Hole for 1/16 " Fitting, PK10
0726632	15% Graphite/85% Vespel Ferrule with 0.2 mm Hole for 1/16 " Fitting, PK10
0726640	15% Graphite/85% Vespel Ferrule with 1 x 0.5 mm Hole & 1 x 0.3 mm Hole for 1/16 " Fitting, PK10
0726641	15% Graphite/85% Vespel Ferrule with 1 x 0.4 mm Hole & 1 x 0.3 mm Hole for 1/16 " Fitting, PK10
0726642	15% Graphite/85% Vespel Ferrule with 1 x 0.5 mm Hole & 1 x 0.4 mm Hole for 1/16 " Fitting, PK10
0726650	15% Graphite/85% Vespel Ferrule with 3 x 0.3 mm Holes for 1/16 " Fitting, PK10
0726660	15% Graphite/85% Vespel Ferrule with 0.5 mm Hole for 1/4 " Fitting, PK10
0726673	15% Graphite/85% Vespel Ferrule with 6 mm Hole for 1/4 " Fitting, PK10
0726700	15% Graphite/85% Vespel Ferrule with 2.4 mm Hole for 1/8 " Fitting, PK10
0726701	15% Graphite/85% Vespel Ferrule with 1/16 " Hole for 1/8 " Fitting, PK10
0726702	15% Graphite/85% Vespel Ferrule with 0.5 mm Hole for 1/8 " Fitting, PK10
0726703	15% Graphite/85% Vespel Ferrule with 0.4 mm Hole for 1/8 " Fitting, PK10
0726288	Replacement Ferrule for Obsolete SGE Varian Flame Jet, PK10
072629	Replacement Mini Union Graphite Ferrules with 0.8 mm Hole, PK10
072630	Replacement Mini Union Graphite Ferrules with 0.5 mm Hole, PK10
072635	Graphite Short Reducing Ferrule with 0.5 mm Hole for Detector /Injectors (Not for MSD Interface Connections), PK10
072636	Graphite Short Reducing Ferrule with 0.8 mm Hole for Detector /Injectors (Not for MSD Interface Connections), PK10

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).