

Tubing | MyCapLC™ Kit

Customize CapLC for Your Application

Part Number	Part Description and Detail
2251110	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 50 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251111	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 50 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251210	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 100 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251211	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 100 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251310	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 150 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251311	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 150 µm ID x 150 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251130	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 50 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251131	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 50 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251230	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 100 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251231	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 100 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251330	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 150 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 360 µm OD x 50 µm ID
2251331	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 300 µm ID x 150 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251151	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 530 µm ID x 50 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251251	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 530 µm ID x 100 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID
2251351	MyCapLC Kit Containing Components to Assemble Three Complete 1/16 " OD x 530 µm ID x 150 mm Length UHPLC Columns, with 250 mm Inlet and 100 mm Outlet Tails 1/32 " OD x 50 µm ID

For more information about this product visit www.sge.com or contact techsupport@sge.com.