

## Syringes | Diamond Headspace

### Higher Limits in Headspace Analysis

Part Number	Part Description and Detail	Replacement Plunger
<b>CTC Head Space</b>		
008131	1 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	031841
008136	1 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.47 mm OD Side-hole Needle	031841
008631	2.5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	031851
008636	2.5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.47 mm OD Side-hole Needle	031851
008831	5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	031861
<b>CTC PAL3/RTC &amp; Thermo RSH Head Space</b>		
008155	1 mL Fixed Needle CTC PAL3/RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 6.5 cm 0.63 mm OD Side-hole Needle	032841
008655	2.5 mL Fixed Needle CTC PAL3/RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 6.5 cm 0.63 mm OD Side-hole Needle	032846
008855	5 mL Fixed Needle CTC PAL3/RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 6.5 cm 0.63 mm OD Side-hole Needle	032851

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).