

Accessories | Hydrocarbon Trap

Protect Your GC Column

Gas purifiers are an essential part of your GC analysis. Without the appropriate gas purifier, problems such as noisy baselines, moisture entering GC column, excessive bleed, reduced column lifetime and irreproducible results will occur. These problems are created by hydrocarbons, moisture and oxygen in the carrier gases.

Features and Benefits


The Hydrocarbon Trap removes C5 weight and heavier hydrocarbons from inert gases, air or hydrogen to low ppb levels. It has an impervious outer stainless steel casing.

Recommended Applications

Recommended for GC carrier or detector gas lines.

Product Specifications

Description	Fitting
Hydrocarbon Trap	1/8 "
Hydrocarbon Trap	1/4 "



For more information about this product visit www.sge.com or contact techsupport@sge.com.