

Tubing | Stainless Steel

SGE stainless steel tubing is available in a number of different size combinations, suitable for your tubing needs.

Features and Benefits

- Smooth and clean bore
- 5, 10 and 25 meter lengths supplied in coils. Where tubing OD is smaller than 1/16", ends are bushed to 1/16" and GC ferrules are supplied.
- 5, 10 and 25 meter lengths supplied cleaned and passivated.

Recommended Applications

GC gas lines.



Product Specifications

OD	ID	Length
1/48"	0.25 mm	25 m
1/32"	0.5 mm	25 m
1/8"	2.2 mm	5 m
1/16"	0.8 mm	10 m
1/16"	0.8 mm	5 m
1/8"	2.2 mm	2 m
1/16"	0.8 mm	2 m

For more information about this product visit www.sge.com or contact techsupport@sge.com.