

Accessories | Syringe Valves

Extraction to Injection in a Single Process

Syringe valves are ideal for sample storage and transportation, or to pressurize sample prior to injection. The valve allows fluid to be stored in the syringe without contamination.

SGE Luer Tip, Luer Lock and removable needle syringes can be fitted with valves.

Available in two configurations of manually operated syringe valves; push-pull for smaller volumes (25 µL – 2.5 mL), and push-button for larger volumes (5 mL – 2 L).

Recommended Applications

Sample storage and transportation. Pressurizing samples prior to injection.



Push-pull Valve with Luer Lock



Push-button Valve with Luer Lock or removable needle



Push-button Valve for Luer Lock devices

Product Specifications

- Push-pull valve with replaceable cone tip needle available for 25 µL - 2.5 mL removable needle syringes.
- Push-button valve with Luer Lock termination available for 5 mL - 2 L syringes with removable Luer Lock or removable needle.
- Push-button Luer Lock valve available for any Luer Lock or Luer Tip syringe, or other Luer devices.

For more information about this product visit www.sge.com or contact techsupport@sge.com.