

Syringes | Guided plunger

Guided plunger syringes are the most rugged syringe available, making them perfect for industrial environments.

Features and benefits

- Extended barrel guides plunger during injection.
- Robust and rugged.

Recommended applications

Ideal for student and industrial use.

Product specifications

Available in 5 µL to 10 µL capacities.

Accuracy and reproducibility	±1% (dispensed volume)
Temperature range	5-70°C fixed needle 5-120°C removable needle



For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com