

Separations | GC Column BPX90

Part Number	Part Description and Detail
054570	0.25 mm ID x 0.25 µm Film Thickness x 15 m Length BPX90 GC Capillary Column
054580	0.25 mm ID x 0.25 µm Film Thickness x 30 m Length BPX90 GC Capillary Column
054590	0.25 mm ID x 0.25 µm Film Thickness x 60 m Length BPX90 GC Capillary Column
054596	0.25 mm ID x 0.25 µm Film Thickness x 100 m Length BPX90 GC Capillary Column
054583	0.32 mm ID x 0.5 µm Film Thickness x 30 m Length BPX90 GC Capillary Column

For more information about this product visit www.sge.com or contact techsupport@sge.com.