

Separations | HPLC Guard Column ProteCol®

Turn Imagination into Results

ProteCol® guard columns are recommended to protect the analytical column, and ensure it performs consistently.

Features and Benefits

The unique design of SGE's guard column couples directly to the analytical column. The guard column is designed to fit into the back of a PEEK fingertight fitting (provided with the guard column). No further unions are required.

Recommended Applications

We recommend the use of ProteCol Guard columns to protect the analytical column, and ensure it performs consistently.

Product Specifications

Available in C18, C8, C4, Phenyl Hexyl, Silica, Amino, Cyano, HILIC, Chiral, SCX and PFP phases.

For more information about this product visit www.sge.com or contact techsupport@sge.com.

