

Injection Port | Inlet Liner Taper FocusLiner®

Spot the Difference in the SGE Lineup!

Part Number	Part Description and Detail
Agilent	
09200301	Agilent Inlet Liner 4 mm ID Tapered FocusLiner, Length 78.5 mm PK1
092003	Agilent Inlet Liners 4 mm ID Tapered FocusLiner, Length 78.5 mm PK5
092011	Agilent Inlet Liners 4 mm ID Tapered FocusLiner, Length 78.5 mm PK25
09211101	Agilent Inlet Liner 2.3 mm ID Tapered FocusLiner, Length 78.5 mm PK1
092111	Agilent Inlet Liners 2.3 mm ID Tapered FocusLiner, Length 78.5 mm PK5
092115	Agilent Inlet Liners 2.3 mm ID Tapered FocusLiner, Length 78.5 mm PK25
Bruker/Varian	
09203601	Bruker/Varian 1075/1077 Inlet Liner 3.4 mm ID Tapered FocusLiner, Length 54 mm PK1
092036	Bruker/Varian 1075/1077 Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 54 mm PK5
09202501	Bruker/Varian 1075/1077 Inlet Liner 4 mm ID Tapered FocusLiner, Length 72 mm PK1
Dani	
092506	Dani Inlet Liners 4 mm ID Tapered FocusLiner, Length 70 mm PK5
PerkinElmer	
09209501	PerkinElmer Inlet Liner 4 mm ID Tapered FocusLiner, Length 92 mm PK1
092095	PerkinElmer Inlet Liners 4 mm ID Tapered FocusLiner, Length 92 mm PK5
09209525	PerkinElmer Inlet Liners 4 mm ID Tapered FocusLiner, Length 92 mm PK25
Shimadzu	
09206801	Shimadzu 2010/17A Inlet Liner 3.4 mm ID Tapered FocusLiner, Length 95 mm PK1
092068	Shimadzu 2010/17A Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 95 mm PK5
09206825	Shimadzu 2010/17A Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 95 mm PK25
09205801	Shimadzu 2010/17A Inlet Liner 3.4 mm ID Tapered FocusLiner, Length 95 mm PK1
092058	Shimadzu 2010/17A Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 95 mm PK5
09205825	Shimadzu 2010/17A Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 95 mm PK25
09206601	Shimadzu 14A Inlet Liner 3.4 mm ID Tapered FocusLiner, Length 99 mm PK1
092066	Shimadzu 14A Inlet Liners 3.4 mm ID Tapered FocusLiner, Length 99 mm PK5

Note - Bruker/Varian 1077 injectors and Thermo Trace 1300 series use 78.5 mm "Agilent" Liners

For more information about this product visit www.sge.com or contact techsupport@sge.com.