

Separations | HPLC Column ProteCol® Phenyl Hexyl

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and Benefits

ProteCol® Phenyl-Hexyl uses a hexyl-linked phenyl phase where the hexyl alkyl chain delivers unique selectivity and increased hydrolytic stability when compared to propyl-linked chemistry.

Recommended Applications

Recommended for compounds that are difficult to separate on standard phenyl phase columns.

Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	pH	USP Code
Phenyl Hexyl	120	3, 5	350	10.0	1-10	L11

For more information about this product visit www.sge.com or contact techsupport@sge.com.

