

Separations | HPLC Column ProteCol® SCX

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and Benefits

- High loading capacity and pressure limit.
- Robust bonding technology.
- High density bonding.
- Bonding aromatic sulfonic acid groups.

Recommended Applications

Small organic bases.

Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

PEEK coated column hardware available to reduce non-specific interactions when analyzing compounds with metal chelating activity.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	pH
SCX	120	3, 5	350	1-10

For more information about this product visit www.sge.com or contact techsupport@sge.com.

