

Accessories | Oxygen Trap

Protect Your GC Column

Gas purifiers are an essential part of your GC analysis. Without the appropriate gas purifier, problems such as noisy baselines, moisture entering GC column, excessive bleed, reduced column lifetime and irreproducible results will occur. These problems are created by hydrocarbons, moisture and oxygen in the carrier gases.

Features and Benefits

- Huge 2.5 L capacity; reduces O₂ to low ppb levels.
- Advanced adsorbent technology.
- Outstanding efficiency in a small, high performance design.

Recommended Applications

Recommended for GC carrier or detector gas lines.

Product Specifications

Description	Fitting
Oxygen Trap	1/8 "
Oxygen Trap	1/4 "

For more information about this product visit www.sge.com or contact techsupport@sge.com.

