


## Connections | HPLC Hexnut™

### Quick and Easy HPLC Connections

Ideal for applications where corrosive solvents are being used, Hexnuts™ have inert contact surfaces, making them biocompatible. Stainless steel 10-32 thread fittings use a non-swaging Kel-F® or PEEK replaceable ferrule.



#### Features and Benefits

- Stainless steel 10-32 thread fittings use a non-swaging Kel-F or PEEK replaceable ferrule.
- Inert contact surfaces make them bio compatible.
- Ideal for applications where corrosive solvents are being used.

#### Recommended Applications

- Biotechnology.
- Medical research.
- Proteomics.
- Pharmaceuticals.
- Recommended for all HPLC connections, particularly where PEEK or PEEKsil® tubing is used.

---

#### Product Specifications

- Maximum Pressure Rating: 5,000 psi
- Thread Type: 10-32

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).