

Accessories | Moisture Trap

Protect Your GC Column

Gas purifiers are an essential part of your GC analysis. Without the appropriate gas purifier, problems such as noisy baselines, moisture entering GC column, excessive bleed, reduced column lifetime and irreproducible results will occur. These problems are created by hydrocarbons, moisture and oxygen in the carrier gases.

Features and Benefits

The Moisture Trap removes water from inert gases, air or hydrogen to less than 20 ppb at flow rates up to 2 L/min.

Recommended Applications

Recommended for GC carrier or detector gas lines.

Product Specifications

The trap has a capacity for water of 15 g.

Description	Fitting
Moisture Trap	1/8 "
Moisture Trap	1/4 "



For more information about this product visit www.sge.com or contact techsupport@sge.com.