

Separations | HPLC Column ProteCol® C4

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and Benefits

High durability and extended acidic and alkaline resistance.

Recommended Applications

Recommended for compounds too strongly retained on C18 and C8 phases, and protein purification.

Product Specifications

0.15, 0.3, 1, 2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths. PEEK coated, and PEEKsil®, column hardware available to reduce non-specific interactions when analyzing compounds with metal chelating activity.

Phase	Pore Size (Å)	Particle Size (µm)	Particle Volume	Surface Area	% Carbon Load	pH	USP Code
C4	300	5	0.9	100	3	1-11	L26

For more information about this product visit www.sge.com or contact techsupport@sge.com.

