

Connections | GC Ferrules Graphite

Part Number	Part Description and Detail
0726001	Graphite Sealing Ring with 5 mm Hole for Shimadzu 14, 15A & 16 (Spl-14 Injector) Liners & 5 mm OD Packed Columns, PK10
072601	Graphite Short Ferrule with 1/4 " Hole for 1/4 " Fittings, Bruker/Varian 1075 & 1077 Injection Port Liners, Unions, & 1/4 " OD Packed Columns, PK10
072602	Graphite Short Ferrule with 1/8 " Hole for 1/8 " Fittings, Unions, & 1/8 " OD Packed Columns, PK10
072603	Graphite Short Ferrule with 1/16 " Hole for 1/16 " Fittings, Unions, PK10
072604	Graphite Short Reducing Ferrule with 6 mm Hole for 1/4 " Fittings, PK10
072605	Graphite Short Ferrule with 5 mm Hole for 5 mm Fittings, PK10
072606	Graphite Short Reducing Ferrule with 1.2 mm Hole for 1/16 " Fittings, PK10
072607	Graphite Short Reducing Ferrule with 1 mm Hole for 1/16 " Fittings, PK10
072608	Graphite Short Reducing Ferrule with 0.8 mm Hole for 1/16 " Fittings, PK10
0726080	Graphite Ferrule with 0.5 mm Hole for GC14A, GC17A, GC2010 and GC2014 Detector / Injectors (Not for MS Interfaces or QP2010 Injector), PK10
0726082	Graphite Ferrule with 0.8 mm Hole for GC14A, GC17A, GC2010 and GC2014 Detector / Injectors (Not for MS Interfaces or QP2010 Injector), PK10
072609	Graphite Short Reducing Ferrule with 1.2 mm Hole for 1/8 " Fittings, PK10
072610	Graphite Short Reducing Ferrule with 0.5 mm Hole for 1/16 " Fittings, PK10
072613	Graphite Short Ferrule with 1/2 " Hole for 1/2 " Fittings, PK10
072614	Graphite Ferrule with 0.8 mm Hole for Carlo Erba Detector / Injectors (Not for GC/MS, Not Suitable for 8000 Series and Trace/Focus GC), PK10
072615	Graphite Short Ferrule with 3/8 " Hole for 3/8 " Fittings, PK10
0726151	Graphite Short Reducing Ferrule with 1/8 " Hole for 1/4 " Fittings, PK10
072617	Graphite Short Reducing Ferrule with 5/16 " Hole for 1/4 " Fittings, PK10
072619	Graphite Ferrule with 0.5 mm Hole for Carlo Erba Detector / Injectors (Not for GC/MS, Not Suitable for 8000 Series and Trace/Focus GC), PK10
072621	Graphite Ferrule with 1/4 " Hole for 1/4 " OD Packed Columns, PK10
0726212	Graphite Reducing Ferrule with 6 mm Hole for 1/4 " Fittings, PK10
0726217	Graphite Sealing Ring for Varian 1078 & 1079 Injectors, PK10
072622	Graphite Ferrule with 1/8 " Hole for 1/8 " OD Packed Columns, PK10
072623	Graphite Ferrule with 1/16 " Hole for 1/16 " Fittings, PK10
072624	Graphite Ferrules with 0.5 mm Hole for Use with 1/8 " Nuts on PerkinElmer Injectors & Detectors, PK10
072625	Graphite Ferrule with 1.0 mm Hole for 1/16 " Fittings, PK10
072626	Graphite Ferrule with 0.8 mm Hole for 1/16 " Fittings, PK10
072627	Graphite Ferrule with 0.5 mm Hole for 1/16 " Fittings, PK10
0726270	Graphite Ferrule with 0.8 mm Hole for 1/8 " Fittings, PK10
072628	Graphite Ferrule with 1.2 mm Hole for 1/8 " Fittings, PK10
0726280	Graphite Ferrules with 1.0 mm Hole for Use with 1/8 " Nuts on PerkinElmer Injectors & Detectors, PK10
0726281	Graphite Ferrule with 1/8 " Hole for 1/4 " Fittings, PK10
0726282	Graphite Ferrule with 1/16 " Hole for 1/4 " Fittings, PK10
0726283	Graphite Ferrule with 1/16 " Hole for 1/8 " Fittings, PK10
0726288	Replacement Ferrule for Obsolete SGE Varian Flame Jet, PK10
072629	Replacement Mini Union Graphite Ferrules with 0.8 mm Hole, PK10
072630	Replacement Mini Union Graphite Ferrules with 0.5 mm Hole, PK10
072635	Graphite Short Reducing Ferrule with 0.5 mm Hole for Detector / Injectors (Not for MSD Interface Connections) , PK10
072636	Graphite Short Reducing Ferrule with 0.8 mm Hole for Detector / Injectors (Not for MSD Interface Connections) , PK10

For more information about this product visit www.sge.com or contact techsupport@sge.com.