

Separations | HPLC Capillary Column ProteCol®

Turn Imagination into Results

The ProteCol® capillary HPLC design includes integrated connection tubing with PEEKsil® for the column body and connection capillaries delivering robust design, zero volume connections and uninterrupted flow.


Features and Benefits

Integrated PEEKsil connection tubing, PEEKsil column body, and zero volume connections allows for use with:

- Small samples.
- Exotic solvents.
- Low concentrations.
- Direct coupling into MS.

Recommended Applications

- Biotechnology.
- Medical research.
- Proteomics.
- Pharmaceuticals.



Product Specifications

150 µm and 300 µm. Available in C18, C8 and C4 phases.

For more information about this product visit www.sge.com or contact techsupport@sge.com.