

Separations | HPLC Column ProteCol® Silica

Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

Features and Benefits

- High durability and extended acidic and alkaline resistance.
- High surface area and mechanical strength.

Recommended Applications

Separation of very hydrophobic molecules.

Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	USP Code
Silica	120	3, 5, 10	300	0	L3

For more information about this product visit www.sge.com or contact techsupport@sge.com.

