

Syringes | Cross Reference Hamilton

Variation Notes:
 FN = Fixed Needle, GT = Gas Tight, with PTFE Tipped plunger, H = Hamilton, LL = Luer Lock, LT = Luer Tip, RN = Removable Needle, WO = Without

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
0160310	1750 HITACHI SYRINGE	007660	500 µL Hitachi LC Autosampler Syringe with GT Plunger with M10 - 1.0 Screw Thread	E	
202640	701RSN AOC 14/17(22s/1.69 "/2)	002898	10 µL Removable Needle Shimadzu Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	23 g Needle
202660	1002LTN CTC SYR (22/56/5)	008631	2.5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	O	23 g Needle
202880	701RN VARIAN 8100/8200	002924	10 µL Fixed Needle CTC/Thermo Syringe with GT Plunger & 5.3 cm 0.5 mm OD Needle	O	25 g Needle
203043	1702N CTC SYRINGE (26s/51/AS)	003700	25 µL Fixed Needle CTC Syringe with GT Plunger & 5 cm 0.47 mm OD Cone Tipped Needle	E	
203072	SYR,701 CTC,SLIM (26s/51/2)	002980	10 µL Fixed Needle CTC/Thermo Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	O	Cone Tipped Needle
203073	701N CTC SYRINGE (22s/51/3)	002710	10 µL Fixed Needle CTC Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
203074	SYR,1702CTC,SLIM (26s/51/AS)	0039871	25 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with GT Plunger & 5 cm 0.63 mm OD Cone Tipped Needle	E	23 g Needle
203075	1702N CTC SYRINGE (22s/51/3)	003715	25 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle (8 mm OD barrel)	E	
203076	SYR,1710 CTC,SLIM (26s/51/AS)	005700	100 µL Fixed Needle CTC Syringe with GT Plunger & 5 cm 0.47 mm OD Cone Tipped Needle	E	
203077	SYR,1710N CTC,SLIM (22s/51/3)	005715	100 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
203078	1725N CTC SYRINGE (26/51/AS)	006700	250 µL Fixed Needle CTC Syringe with GT Plunger & 5 cm 0.47 mm OD Cone Tipped Needle	E	
203082	1001LTN CTC SYRINGE (23/51/5)	008131	1 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	E	
203084	1002LTN CTC SYRINGE (23/51/5)	008631	2.5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	E	
203141	1001LTN CTC SYRINGE (26/51/5)	008136	1 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.47 mm OD Side-hole Needle	E	
203181	1002LTN CTC SYRINGE (26/51/5)	008636	2.5 mL Fixed Needle CTC Head Space Syringe with Energized GT Plunger and 5.6 cm 0.47 mm OD Side-hole Needle	E	
203185	7701.2 CTC SYRINGE (26/51/AS)	000790	2.0 µL NanoVolume CTC/Thermo Syringe with 5.0 cm 0.63 mm OD Cone Tipped Needle	O	2 µL Instead of 1 µL
203189	75N CTC SYRINGE (26s/51/AS)	001982	5 µL Fixed Needle CTC/Thermo Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
203205	701N CTC SYRINGE (26s/51/AS)	002700	10 µL Fixed Needle CTC Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
203235	1710N CTC SYRINGE (22/51/3)	005720	100 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD (0.4 mm ID) LC Needle	E	
203349	1750N CTC SYR (22/51/3)	007720	500 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD (0.4 mm ID) LC Needle	E	
203361	701N CTC SYR (23S/51/AS)	002981	10 µL Fixed Needle CTC/Thermo Syringe with 5 cm 0.63 mm OD Cone Tipped Needle	E	
203362	701N CTC SYR (23S-26S/51/AS)	002980	10 µL Fixed Needle CTC/Thermo Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	O	26 g Needle
203363	701N CTC SYR (23S/2 "/2)	002981	10 µL Fixed Needle CTC/Thermo Syringe with 5 cm 0.63 mm OD Cone Tipped Needle	O	Cone Tipped Needle
203566	1710 CTC A200S SYR (22s/2/3)	005331	100 µL Fixed Needle CTC/Thermo Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
204000	SYRINGE,75N(26S/50/2)	0019821	5 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	O	Cone Tipped Needle
204001	SYRINGE,701SN(26S/50/2)	0029801	10 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	O	Cone Tipped Needle
204051	SYRINGE,75SN(26S/50/AS)	0019821	5 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
204052	SYRINGE,701SN 26S/50/AS)	0029801	10 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
207189	701 PAL3 (565/85/AS)	002867	10 µL Fixed Needle CTC RTC & Thermo RSH Syringe with 8.5 cm 0.47 mm OD Cone Tipped Needle	E	
207290	1002 CTC PAL3 (23/65/5)HS	008655	2.5 mL Fixed Needle CTC RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 5.6 cm 0.63 mm OD Side-hole Needle	E	
207292	1010 CTC PAL3(19/57/3)	008951	10 mL Removable Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 1.07 mm OD LC Tipped Needle	E	
207803	701 PAL3 (235/57/AS)	002861	10 µL Fixed Needle CTC RTC & Thermo RSH Syringe with 5.7 cm 0.63 mm OD Cone Tipped Needle	E	
207807	701 PAL3 (265/57/AS)	002865	10 µL Fixed Needle CTC RTC & Thermo RSH Syringe with 5.7 cm 0.47 mm OD Cone Tipped Needle	E	
207813	1701 PAL3 (225/57/3)	002880	10 µL Removable Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.72 mm OD LC Needle	E	
207829	1725 PAL3 (23/57/AS)	006862	250 µL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.63 mm OD Cone Tipped Needle	E	
207831	1725 PAL3 (26/85/AS)	006862	250 µL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.63 mm OD Cone Tipped Needle	E	
207836	1001 PAL3 (22/57/3)	008150	1 mL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.72 mm OD LC Tipped Needle	E	
207837	1701 PAL3 (265/57/AS)	002866	10 µL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.47 mm OD Cone Tipped Needle	E	
207839	1701 PAL3 (235/57/AS)	002862	10 µL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.63 mm OD Cone Tipped Needle	E	
207840	1710 PAL3 (265/57/AS)	005866	100 µL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 0.47 mm OD Cone Tipped Needle	E	
207853	1001 CTC PAL3 (23/65/5)HS	008155	1 mL Fixed Needle CTC RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 6.5 cm 0.63 mm OD Side-hole Needle	E	
208311	1710 CTC SYR X PAL3 (22/57/3)	008950	10 mL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 1.07 mm OD LC Tipped Needle	E	
208312	1710 CTC SYR X PAL3 (225/57/3)	008950	10 mL Fixed Needle CTC RTC & Thermo RSH Syringe with GT Plunger & 5.7 cm 1.07 mm OD LC Tipped Needle	E	
54658-01	SYRINGE,1710C,AUTO SAMPLER	005696	100 µL V6 Dispenser Syringe	E	
54659-01	SYRINGE,1725C,AUTO SAMPLER	006696	250 µL V6 Dispenser Syringe	E	
54660-01	SYRINGE, 1750C,AUTO SAMPLER	007696	500 µL V6 Dispenser Syringe	E	
54661-01	SYRINGE,1001C,AUTO SAMPLER	008196	1.0 mL V6 Dispenser Syringe	E	
54662-01	SYRINGE,1002C,AUTO SAMPLER	008596	2.5 mL V6 Dispenser Syringe	E	
62161-01	SYRINGE,1702C,SPARK AUTO	003696	25 µL V6 Dispenser Syringe	E	
67430-01	SYR,1702N,(23/51/AS),CTC-S	003700	25 µL Fixed Needle CTC Syringe with GT Plunger & 5 cm 0.47 mm OD Cone Tipped Needle	E	
67432-01	SYR,1710N,(225/51/3),CTC-S-X	005331	100 µL Fixed Needle CTC/Thermo Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
67434-01	SYR,1710N,(23/51/AS),CTC-S	005335	100 µL Fixed Needle CTC/Thermo Syringe with GT Plunger & 5 cm 0.63 mm OD Cone Tipped Needle	E	
67436-01	SYR,75N,(265/51/AS),CTC-S	0019821	5 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
67438-01	SYR,701N,(265/51/AS),CTC-S	0029801	10 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
67440-01	SYR,701SN,(235/51/AS),CTC-S	002951	10 µL Fixed Needle Varian 8400/CTC/Thermo (Classic Button) Syringe with 5 cm 0.63 mm OD Cone Tipped Needle	E	
67442-01	SYR,1725N,(22/51/3),CTC-S	006330	250 µL Removeable Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD (0.4 mm ID) LC Needle	E	
67446-01	SYR,1702N,(225/51/3),CTC-S	003715	25 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle (8 mm OD barrel)	E	
67448-01	SYR,1750N,(22/51/3),CTC-S	007720	500 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD (0.4 mm ID) LC Needle	E	
67450-01	SYR,1705N,(22/51/3),CTC-S	004810	50 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
67452-01	SYR,1710SN,(225/51/3),CTC-S	005715	100 µL Fixed Needle CTC Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
67454-01	SYR,1701SN,(235/51/AS),CTC-S	0029871	10 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with GT Plunger & 5 cm 0.63 mm OD Cone Tipped Needle	E	
80000	1701N 10 µL SYR (26s/2 "1/2)	002200	10 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.47 mm OD Bevel Tipped Needle	E	
80004	1701NCH 10 µL SYR (26s/2 "1/2)	002200 + 031930	10 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.47 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80020	1702 WISP SYR 25 µL	003990	25 µL Waters Wisp Syringe with GT Plunger	E	
80024	1725 WISP SYR 250 µL	006690	250 µL Waters Wisp Syringe with GT Plunger	E	
80030	1701RN 10 µL SYR (26s/2 "1/2)	002250	10 µL Removable Needle Syringe with GT Plunger & 5 cm 0.47 mm OD Bevel Tipped Needle	E	
80039	1701NPT5 10 µL SYR (26s/2 "1/5)	002250 + 037410	10 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
80065	1701RNR 10 µL SYR (22s/2 °/3)	002313	10 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
80074	SYRINGE,175ASN,23S GA,1.71	001810	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	Not GT (Note that H is Half Scale 10 µL)
80074	SYRINGE,175ASN,23S GA,1.71	002812	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63 mm OD Cone Tipped Needle	O	Full Scale 10 µL Instead of Half Scale 5 µL
80075	1701N 10 µL SYR (26s/2 °/3)	002335	10 µL Fixed Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	22 g Needle
80076	175ASN SYR (23-26/1.71 °/HP)	001821	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Not GT (Note that H is Half Scale 10 µL)
80076	175ASN SYR (23-26/1.71 °/HP)	002826	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Full Scale 10 µL Instead of Half Scale 5 µL
80079	1701ASN SYR(23S/26S/1.71 °/HP)	002826	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	E	
80080	1701ASN SYR (23s/1.71 °/HP)	002812	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
80084	175ASRN(23s/1.71 °/HP)	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Dual Gauge Needle (Note that H 5 µL is Half Scale 10 µL)
80085	1701N 10 µL SYR (26/2 °/3)	002335	10 µL Fixed Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	22 g Needle
80086	175ASRN SYR (23-26/1.71 °/HP)	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Full Scale 10 µL Instead of Half Scale 5 µL
80087	1701ASRN 10 µL SYR 23s/1.71 °/HP	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Dual Gauge Needle
80088	1701ASRN 10 µL SYR 26s/1.71 °/HP	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	Dual Gauge Needle
80089	1701ASRN SYR (23S/26S/1.71/HP)	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	E	
80090	175ASN 5 µL (23S/1.71 °/HP) 6/PK	002812	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63 mm OD Cone Tipped Needle	O	Full Scale 10 µL Instead of Half Scale 5 µL
80090	175ASN 5 µL (23S/1.71 °/HP) 6/PK	002816	10 µL Removable Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle PK25	O	Full Scale 10 µL Instead of Half Scale 5 µL
80092	175ASN 5 µL (23-26/1.71 °/HP)6/PK	002827	5 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle PK6	O	Full Scale 10 µL Instead of Half Scale 5 µL
80094	1701ASN SYR (23s/1.71 °/HP)6/PK	002812	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63 mm OD Cone Tipped Needle	O	
80094	1701ASN SYR (23s/1.71 °/HP)6/PK	002816	10 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63 mm OD Cone Tipped Needle PK25	O	
80096	1701ASN (23-26/1.71 °/HP) 6/PK	002827	5 µL Fixed Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle PK6	E	
80100	7001KH 1 µL SYR (25/2.75 °/3)	000570	1.0 µL NanoVolume Syringe with 7 cm 0.47 mm OD Cone Tipped Needle	O	26 g, Cone Tipped Needle
80104	7001KHCH 1 µL SYR (25/2.75 °/3)	000570 + 031930	1.0 µL NanoVolume Syringe with 7 cm 0.47 mm OD Cone Tipped Needle and Repeating Adaptor	O	26 g, Cone Tipped Needle
80107	7001KHWG 1 µL SYR (25s/2.75 °/3)	000570 + 031930	1.0 µL NanoVolume Syringe with 7 cm 0.47 mm OD Cone Tipped Needle and Repeating Adaptor	O	26 g, Cone Tipped Needle
80135	7001KH 1.0 µL SYR (25/2.75 °/2)	000506	1.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Bevel Tipped Needle	O	23 g Needle
80175	7001ASRN (26/1.71 °/HP)	000610	1.0 µL NanoVolume Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	23 g Needle
80176	7001ASRN (23/1.71 °/HP)	000610	1.0 µL NanoVolume Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
80200	1702N 25 µL SYR (22s/2 °/2)	003200	25 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	O	25 g Needle
80201	1702LT 25 µL SYR	004229	50 µL Fixed Luer Tip Syringe with GT Plunger	O	50 µL
80222	1702TLLX 25 µL SYR	004232	50 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	O	50 µL
80223	1702TLLX SYR W/BSHG STD BTN	004232	50 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	O	50 µL
80230	1702RN 25 µL SYR (22s/2 °/2)	003250	25 µL Removable Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	O	25 g Needle
80238	1702RNW 25 µL SYR(25s/1.97 °/3)	003250 + 038260	25 µL Removable Needle Syringe with GT Plunger & Waters LC Tipped Needle	E	
80239	1702NPT5 25 µL SYR (22s/2 °/5)	003250 + 038410	25 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80265	1702RNR 25 µL SYR (22s/2 °/3)	003312	25 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
80275	1702N 25 µL SYR (22s/2 "/3)	003300	25 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	O	Not GT
80285	1702N 25 µL SYR (22/2 "/3)	003300	25 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	O	Not GT
80300	701N 10 µL SYR (26s/2 "/2)	002000	10 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	Manual Use
80300	701N 10 µL SYR (26s/2 "/2)	002950	10 µL Fixed Needle Varian 8400 Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	Bruker/Varian Autosamplers
80304	701NCH 10 µL SYR (26s/2 "/2)	002000 + 031930	10 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80307	701NWG 10 µL SYR (26s/2 "/2)	002000 + 031930	10 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80318	701SN 10 µL SYR (26s/3.15 "/HP)	0029921	10 µL Fixed Needle CTC/Thermo (Classic Button) Syringe with 8 cm 0.47 mm OD Cone Tipped Needle	E	
80330	701RN 10 µL SYR (26s/2 "/2)	002050	10 µL Removable Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	
80334	701RNCH 10 µL SYR (26s/2 "/2)	002050 + 031930	10 µL Removable Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80336	701RN 10 µL SYR 6/PK (26s/2 "/2)	002080	10 µL Removable Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle PK6	E	
80337	701RNWG 10 µL SYR (26s/2 "/2)	002050 + 031930	10 µL Removable Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80339	701NPT5 10 µL SYR (26s/2 "/5)	002050 + 037410	10 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN
80342	701SN 10 µL SYR (23s/2 "/HP)	002005	10 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Cone Tipped Needle	E	
80350	701SN 10 µL SYR (26s/2.75 "/2)	002003	10 µL Fixed Needle Syringe with 7 cm 0.47 mm OD Bevel Tipped Needle	E	
80357	701ASRN 10 µL SYR (23s/1.71 "/HP)	002815	10 µL Removable Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
80358	701ASRN 10 µL SYR (26s/1.71 "/HP)	002805	10 µL Removable Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	E	
80359	701ASRN 10 µL (23-26/1.71 "/HP)	002829	10 µL Removable Needle Agilent Syringe with GT Plunger & 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	GT
80360	901N 10 µL SYR (26s/2 "/2)	002400	10 µL Fixed Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	O	
80365	701SNR 10 µL SYR (22s/2 "/3)	002301	10 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80366	701N 10 µL SYR 6/PK (26s/2 "/2)	002030	10 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle PK6	E	Manual Use
80370	901RN 10 µL SYR (26s/2 "/2)	002450	10 µL Removable Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	
80383	701SN 10 µL SYR (26s/2 "/3)	002301	10 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80384	701SN 10 µL SYR (26s/3 "/2)	002003	10 µL Fixed Needle Syringe with 7 cm 0.47 mm OD Bevel Tipped Needle	O	7 cm Needle
80387	701ASN 10 µL SYR (23s/1.71 "/HP)	002810	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
80388	701ASN 10 µL SYR (26s/1.71 "/HP)	002800	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	E	
80389	701ASN 10 µL SYR-HP 26s (6/PK)	002804	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle PK6	E	
80390	701ASN 10 µL SYR-HP 23s (6/PK)	002814	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle PK6	E	
80391	701ASN 6/PK (23-26/1.71 "/HP)	002822	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle PK6	E	
80393	701ASN SYR (23-26/1.71 "/HP)	002821	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	E	
80398	701ASN 10 µL SYR (23s/1.71 "/2)	002810	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	Cone Tipped Needle
80399	701ASN 10 µL SYR (26s/1.71 "/2)	002800	10 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	O	Cone Tipped Needle
80400	702N 25 µL SYR (22s/2 "/2)	003000	25 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80404	702NCH 25 µL SYR (22s/2 "/2)	003000 + 031930	25 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80407	702NWG 25 µL SYR (22s/2 "/2)	003000 + 031930	25 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80419	702SN 25 µL SYR (22/2 "/3)	003300	25 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80430	702RN 25 µL SYR (22s/2 "/2)	003050	25 µL Removable Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80439	702NPT5 25 µL SYR (22s/2 "/5)	003050 + 038410	25 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80465	702SNR 25 µL SYR (22s/2 "/3)	003300	25 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80500	705N 50 µL SYR (22s/2 "/2)	004000	50 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
80501	705LT 50 µL SYR	004229	50 µL Fixed Luer Tip Syringe with GT Plunger	E	
80505	705NCH 50 µL SYR (22s/2 °/2)	004000 + 031930	50 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80521	705SN 50 µL SYR (22/2 °/3)	004300	50 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80530	705RN 50 µL SYR (22s/2 °/2)	004050	50 µL Removable Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80539	705NPT5 50 µL SYR (22s/2 °/5)	004050 + 038410	50 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80565	705SNR 50 µL SYR (22s/2 °/3)	004300	50 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80600	710N 100 µL SYR (22s/2 °/2)	005000	100 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80601	710LT 100 µL SYR	005229	100 µL Fixed Luer Tip Syringe with GT Plunger	E	
80605	710NCH 100 µL SYR (22s/2 °/2)	005000 + 031930	100 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle and Repeating Adaptor	E	
80621	710SN 100 µL SYR (22/2 °/3)	005300	100 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80630	710RN 100 µL SYR (22s/2 °/2)	005050	100 µL Removable Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80639	710NPT5 100 µL SYR (22s/2 °/5)	005050 + 038410	100 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80665	710SNR 100 µL SYR (22s/2 °/3)	005300	100 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80700	725N 250 µL SYR (22S/2 °/2)	006000	250 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80701	725LT 250 µL SYR	006229	250 µL Fixed Luer Tip Syringe with GT Plunger	E	
80730	725RN 250 µL SYR (22s/2 °/2)	006050	250 µL Removable Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80739	725N 250 µL SYR (22S/2 °/5)	006050 + 038410	250 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80765	725NR 250 µL (22/2 °/3)	006300	250 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
80800	750N 500 µL SYR (22/2 °/2)	007000	500 µL Fixed Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80801	750LT 500 µL SYR	007229	500 µL Fixed Luer Tip Syringe with GT Plunger	E	
80830	750RN 500 µL SYR (22/2 °/2)	007050	500 µL Removable Needle Syringe with 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80839	750N 500 µL SYR (22/2 °/5)	007050 + 038410	500 µL Removable Needle Syringe with 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80865	750N 500 µL SYR (22/2 °/3)	007300	500 µL Fixed Needle Syringe with & 5.1 cm 0.72 mm OD LC Needle	E	
80900	1705N 50 µL SYR (22s/2 °/2)	004200	50 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80901	1705LT 50 µL SYR	004229	50 µL Fixed Luer Tip Syringe with GT Plunger	E	
80920	1705TLL 50 µL SYR	004230	50 µL Fixed Luer Lock Syringe with GT Plunger	E	
80922	1705TLLX 50 µL SYR	004232	50 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	E	
80927	1705TLL 50 µL SYR W/SLOTS	004230	50 µL Fixed Luer Lock Syringe with GT Plunger	E	
80930	1705RN 50 µL SYR (22s/2 °/2)	004250	50 µL Removable Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	
80938	1705RNW 50 µL SYR(25s/1.97 °/3)	004250 + 038260	50 µL Removable Needle Syringe with GT Plunger & Waters LC Tipped Needle	E	
80939	1705NPT5 50 µL SYR (22s/2 °/5)	004250 + 038410	50 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle
80956	1705SL 50 µL SYR (22s/2 °/2)L	004279	50 µL Syringe with Removable Needle Valve & GT Plunger	E	
80962	1705CX 50 µL SYR W/STOP	004995	50 µL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
80965	1705RNR 50 µL SYR (22s/2 °/3)	004312	50 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
80975	1705N 50 µL SYR (22s/2 °/3)	004312	50 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
80985	1705N 50 µL SYR (22/2 °/3)	004312	50 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81000	1710N 100 µL SYR (22s/2 °/2)	005200	100 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	O	25 g Needle
81001	1710LT 100 µL SYR	005229	100 µL Fixed Luer Tip Syringe with GT Plunger	E	
81020	1710TLL 100 µL SYR	005230	100 µL Fixed Luer Lock Syringe with GT Plunger	E	
81022	1710TLLX 100 µL SYR	005232	100 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	E	
81027	1710TLL 100 µL SYR W/SLOTS	005230	100 µL Fixed Luer Lock Syringe with GT Plunger	E	
81030	1710RN 100 µL SYR (22s/2 °/2)	005250	100 µL Removable Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	
81038	1710RNW 100 µL SYR(25s/1.97 °/3)	005250 + 038260	100 µL Removable Needle Syringe with GT Plunger & Waters LC Tipped Needle	E	
81039	1710NPT5 100 µL SYR (22s/2 °/5)	005250 + 038410	100 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	RN, 25 g Needle

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
81056	1710SL 100 µL SYR (22s/2 *1/2)L	005279	100 µL Syringe with Removable Needle Valve & GT Plunger	E	
81060	1710C 100 µL SYR	005990	100 µL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81062	1710CX 100 µL SYR W/STOP	005990	100 µL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81065	1710RNR 100 µL SYR (22s/2 *1/3)	005312	100 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
81075	1710N 100 µL SYR (22s/2 *1/3)	005312	100 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81085	1710N 100 µL SYR (22/2 *1/3)	005312	100 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81100	1725LTN 250 µL SYR (22s/2 *1/2)	006200	250 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	O	Not LT
81101	1725LT 250 µL SYR	006229	250 µL Fixed Luer Tip Syringe with GT Plunger	E	
81120	1725TLL 250 µL SYR	006230	250 µL Fixed Luer Lock Syringe with GT Plunger	E	
81122	1725TLLX 250 µL SYR	006232	250 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	E	
81127	1725TLL 250 µL SYR W/SLOTS	006230	250 µL Fixed Luer Lock Syringe with GT Plunger	E	
81130	1725RN 250 µL SYR (22s/2 *1/2)	006250	250 µL Removable Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	
81138	1725RNW 250 µL SYR(25s/1.97 *1/3)	006250 + 038260	250 µL Removable Needle Syringe with GT Plunger & Waters LC Tipped Needle	O	RN
81139	1725LTN 250 µL SYR (22S/2 *1/5)	006250 + 038410	250 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	
81156	1725SL 250 µL SYR (22s/2 *1/2)L	006279	250 µL Syringe with Removable Needle Valve & GT Plunger	E	
81162	1725CX 250 µL SYR W/STOP	006995	250 µL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81165	1725RNR 250 µL SYR (22/2 *1/3)	006312	250 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
81175	1725LTN 250 µL SYR (22S/2 *1/3)	006312	250 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81185	1725LTN 250 µL SYR (22/2 *1/3)	006312	250 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81201	1750LT 500 µL SYR	007229	500 µL Fixed Luer Tip Syringe with GT Plunger	E	
81216	1750LTN 500 µL SYR (22/2 *1/3)	007312	500 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81217	1750LTN 500 µL SYR (22/2 *1/2)	007200	500 µL Fixed Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	Not LT
81220	1750TLL 500 µL SYR	007230	500 µL Fixed Luer Lock Syringe with GT Plunger	E	
81222	1750TLLX 500 µL SYR	007232	500 µL Fixed Luer Lock Syringe with GT Plunger & Plunger Stop	E	
81227	1750TLL 500 µL SYR W/SLOTS	007230	500 µL Fixed Luer Lock Syringe with GT Plunger	E	
81230	1750RN 500 µL SYR (22/2 *1/2)	007250	500 µL Removable Needle Syringe with GT Plunger & 5 cm 0.5 mm OD Bevel Tipped Needle	E	
81243	1750LTNPT5 500 µL SYR (22/2 *1/5)	007250 + 038410	500 µL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	Not LT, 25 g Needle
812560	1750SL	007279	500 µL Syringe with Removable Needle Valve & GT Plunger	E	23 g Needle
81262	1750CX 500 µL SYR W/STOP	007995	500 µL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81265	1750RNR 500 µL SYR (22/2 *1/3)	007312	500 µL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
81301	1001LT 1.0 mL SYR	008020	1 mL Fixed Luer Tip Syringe with GT Plunger	E	
81316	1001LTN 1.0 mL SYR (22/2 *1/3)	008105	1 mL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81317	1001LTN 1.0 mL SYR (22/2 *1/2)	008102	1 mL Fixed Needle Syringe with GT Plunger & 5 cm 0.72 mm OD Bevel Tipped Needle	O	Not LT
81318	1001LTNCH 1.0 mL SYR (22/2 *1/2)	008102	1 mL Fixed Needle Syringe with GT Plunger & 5 cm 0.72 mm OD Bevel Tipped Needle	O	WO Adaptor
81320	1001TLL 1.0 mL SYR	008025	1 mL Fixed Luer Lock Syringe with GT Plunger	E	
81324	1001TLLCH 1.0 mL SYR	008025	1 mL Fixed Luer Lock Syringe with GT Plunger	O	WO Adaptor
81327	1001TLL 1.0 mL SYR W/SLOTS	008025	1 mL Fixed Luer Lock Syringe with GT Plunger	E	
81330	1001RN 1.0 mL SYR (22/2 *1/2)	008100	1 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	E	
81343	1001LTNPT5 1.0 mL SYR (22/2 *1/5)	008100 + 039120	1 mL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	Not LT, 23 g Needle
81356	1001SL 1.0 mL SYR (22/2 *1/2)L	008110	1 mL Syringe with Removable Needle Valve & GT Plunger	E	
81360	1001C 1.0 mL SYR	008185	1 mL Fixed Needle CTC RTC & Thermo RSH Head Space Syringe with Energized GT Plunger and 6.5 cm 0.63 mm OD Side-hole Needle	E	
81365	1001RNR 1.0 mL SYR (22/2 *1/3)	008105	1 mL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	E	
81401	1002LT 2.5 mL SYR	008420	1.0 mL V6 Dispenser Syringe	E	
81403	1002LTCH 2.5 mL SYR	008020	1 mL Fixed Luer Tip Syringe with GT Plunger	O	WO Adaptor
81416	1002LTN 2.5 mL SYR (22/2 *1/3)	008505	2.5 mL Removable Needle Syringe with GT Plunger & 5.1 cm 0.72 mm OD LC Needle	O	RN
81417	1002LTN 2.5 mL SYR (22/2 *1/2)	008502	2.5 mL Fixed Needle Syringe with GT Plunger & 5 cm 0.72 mm OD Bevel Tipped Needle	E	
81418	1002LTNCH 2.5 mL SYR (22/2 *1/2)	008502	2.5 mL Fixed Needle Syringe with GT Plunger & 5 cm 0.72 mm OD Bevel Tipped Needle	O	WO Adaptor
81420	1002TLL 2.5 mL SYR	008425	2.5 mL Fixed Luer Tip Syringe with GT Plunger	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
81424	1002TLLCH 2.5 mL SYR	008425	2.5 mL Fixed Luer Lock Syringe with GT Plunger	O	WO Adaptor
81427	1002TLL 2.5 mL SYR W/SLOTS	008425	2.5 mL Fixed Luer Lock Syringe with GT Plunger	E	
81430	1002RN 2.5 mL SYR (22/2 °/2)	008500	2.5 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	E	
81443	1002LTNP5 2.5 mL SYR (22/2 °/5)	008500 + 039120	2.5 mL Removable Needle Syringe with GT Plunger & 5 cm Side Hole Dome Tipped Needle	O	Not LT, 23 g Needle
81456	1002SL 2.5 mL SYR (22/2 °/2)L	008510	2.5 mL Syringe with Removable Needle Valve & GT Plunger	E	
81460	1002C 2.5 mL SYR	008687	2.5 mL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81501	1005LT 5.0 mL SYR	008760	5 mL Removable Luer Lock Syringe with GT Plunger	O	LL Not LT
81516	1005LTN 5 mL SYR (22/2 °/3)	008700 + 0315233	5 mL Removable Needle Syringe with GT Plunger & LC Needle	O	RN
81517	1005LTN 5.0 mL SYR (22/2 °/2)	008700	5 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	O	RN
81518	1005LTNCH 5.0 mL SYR (22/2 °/2)	008762	5 mL Fixed Luer Lock Syringe with GT Plunger	O	LL Not LT. WO Adaptor
81520	1005TLL 5.0 mL SYR	008762	5 mL Fixed Luer Lock Syringe with GT Plunger	E	
81524	1005TLLCH 5.0 mL SYR	008762	5 mL Fixed Luer Lock Syringe with GT Plunger	O	WO Adaptor
81527	1005TLL 5 mL SYR W/SLOTS	008762	5 mL Fixed Luer Lock Syringe with GT Plunger	E	
81530	1005RN 5.0 mL SYR (22/2 °/2)	008700	5 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	E	
81560	1005C 5.0 mL SYR	008787	5 mL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
81601	1010LT 10.0 mL SYR	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	O	LL Not LT
81610	1010W 10.0 mL (WATERS) SYR	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	E	
81616	1010LTN 10 mL SYR (22/2 °/3)	008900 + 0315233	10 mL Removable Needle Syringe with GT Plunger & LC Needle	O	RN
81617	1010LTN 10.0 mL SYR (22/2 °/2)	008900	10 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	O	RN
81618	1010LTNCH 10.0 mL SYR (22/2 °/2)	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	O	LL Not LT. WO Adaptor
81620	1010TLL 10.0 mL SYR	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	E	
81624	1010TLLCH 10.0 mL SYR	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	O	WO Adaptor
81627	1010TLL 10 mL SYR W/SLOTS	008962	10 mL Fixed Luer Lock Syringe with GT Plunger	E	
81630	1010RN 10.0 mL SYR (22/2 °/2)	008900	10 mL Removable Needle Syringe with GT Plunger & 5 cm 0.63 mm OD Bevel Tipped Needle	E	
81660	1010C 10.0 mL SYR	008987	10 mL Syringe with GT Plunger and 1/4-28 UNF Termination	E	
82017	1001.25LTN 1.25 mL SYR (22/2 °/2)	009462	25 mL Removable Luer Lock Syringe with GT Plunger	O	LL. Select Needle Separately
82520	1025TLL 25.0 mL SYR	009463	25 mL Fixed Luer Lock Syringe with GT Plunger	E	
82525	1025TLL W/O FL 25 mL SYR W/ SLOT	009463	25 mL Fixed Luer Lock Syringe with GT Plunger	O	WO Adaptor
82527	1025TLL 25 mL SYR W/SLOTS	009463	25 mL Fixed Luer Lock Syringe with GT Plunger	E	
84850	85N 5 µL SYR (26s/2 °/2)	001400	5 µL Fixed Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	O	
84851	85RN 5 µL SYR (26s/2 °/2)	001450	5 µL Removable Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	O	
84852	801N 10 µL SYR (26s/2 °/2)	002400	10 µL Fixed Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	O	
84853	801RN 10 µL SYR (26s/2 °/2)	002450	10 µL Removable Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	O	
85020	1050TLL 50.0 mL SYR	009660	50 mL Removable Luer Lock Syringe with GT Plunger	E	
85027	1050TLL 50 mL SYR W/SLOTS	009660	50 mL Removable Luer Lock Syringe with GT Plunger	E	
86020	1100TLL 100.0 mL SYR	009760	100 mL Removable Luer Lock Syringe with GT Plunger	E	
86204	7101KHCH 1 µL SYR (22/2.75 °/3)	000553	1.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	O	Cone Tip
86207	7101KHWG 1 µL SYR (22/2.75 °/3)	000553	1.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	O	Cone Tip
86211	7101KH 1 µL SYR (22s/2.75 °/2)	000505	1.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle	O	23 g Needle
86250	7000.5KH 0.5 µL SYR (25/2.75 °/3)	000310	0.5 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle	O	Cone Tip
86252	7000.5KHCH .5 µL SYR(25/2.75/3)	000355	0.5 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	O	Cone Tip
86254	7000.5KHWG .5 µL SYR(25/2.75/3)	000355	0.5 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	O	Cone Tip
86259	7000.5KH 0.5 µL SYR (25/2.75 °/2)	000311	0.5 µL NanoVolume Syringe with 7 cm 0.63 mm OD Bevel Tipped Needle	O	23 g Needle
86274	7000.5ASRN (26/1.71 °/HP)	000400	0.5 µL NanoVolume Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	E	
86276	7000.5ASRN (23/1.71 °/HP)	000410	0.5 µL NanoVolume Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
86311	S0500 0.5L SYR (TLL)	009910	500 mL Jumbo Luer Lock Syringe	E	
86312	S1000 1.0L SYR (TLL)	009920	1000 mL Jumbo Luer Lock Syringe	E	
86313	S1500 1.5L SYR (TLL)	009930	2000 mL Jumbo Luer Lock Syringe	O	2 L Instead of 1.5 L
86314	S2000 2.0L SYR (TLL)	009930	2000 mL Jumbo Luer Lock Syringe	E	

Syringes | Cross Reference Hamilton Continued

Hamilton Part Number	Part Description and Detail	SGE Part Number	Part Description and Detail	Equivalent or Option	Variation
86326	1025SL 25 mL SYR (22/2 "/2)L	009472	25 mL Syringe with Removable Luer Lock Valve & GT Plunger	O	LL. Select Needle Separately
863260	1025SL	009472	25 mL Syringe with Removable Luer Lock Valve & GT Plunger	O	Select Needle Separately
86336	1050SL 50 mL SYR (22/2 "/2)L	009670	50 mL Syringe with Removable Luer Lock Valve & GT Plunger	O	LL. Select Needle Separately
863360	1050SL	009670	50 mL Syringe with Removable Luer Lock Valve & GT Plunger	O	Select Needle Separately
86346	1100SL 100 mL SYR (22/2 "/2)L	009770	100 mL Syringe with Removable Luer Lock Valve & GT Plunger	O	LL. Select Needle Separately
87402	701RNFS 10 µL SYR(.17mm/10cm/3)	002050 + 037610	10 µL Removable Needle Syringe with 10 cm 0.17 mm OD Needle	E	
87404	1701RNFS 10 µL SYR(.17mm/10cm/3)	002250 + 037610	10 µL Removable Needle Syringe with GT Plunger & 5 cm 0.47 mm OD Bevel Tipped Needle	E	
87900	75N 5 µL SYR (26s/2 "/2)	001000	5 µL Fixed Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	
87919	75N 5 µL SYR (26s/2 "/3)	001301	5 µL Fixed Needle Syringe with 5.1 cm 0.72 mm OD LC Needle	E	
87930	75RN 5.0 µL SYR (26s/2 "/2)	001050	5 µL Removable Needle Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	
87943	65RNR 5.0 µL SYR (22s/2 "/3)	001450	5 µL Removable Needle Guided Plunger Syringe with 5 cm 0.47 mm OD Bevel Tipped Needle	E	
87957	75ASRN 5 µL SYR (23s/1.71 "/HP)	001810	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	FN
87958	75ASRN 5 µL SYR (26s/1.71 "/HP)	001800	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	O	FN
87959	75ASRN SYR (23S/26S/1.71 "/HP)	001821	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	O	FN
87987	75ASN 5 µL SYR (23s/1.71 "/HP)	001810	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	E	
87988	75ASN 5 µL SYR (26s/1.71 "/HP)	001800	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	E	
87989	75ASN 5 µL SYR(26s/1.71 "/HP)6PK	001804	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle PK6	E	
87990	75ASN 5 µL SYR(23s/1.71 "/HP)6PK	001814	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle PK6	E	
87991	75ASN 5 µL SYR (23s/1.71 "/2)	001810	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63 mm OD Cone Tipped Needle	O	Cone Tipped Needle
87992	75ASN 5 µL SYR (26s/1.71 "/2)	001800	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.47 mm OD Cone Tipped Needle	O	Cone Tipped Needle
87993	75ASN 5 µL SYR (23-26/1.71 "/HP)	001821	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle	E	
87994	75ASN 5 µL 6/PK(23-26/1.71 "/HP)	001822	5 µL Fixed Needle Agilent Syringe with 4.2 cm 0.63/0.47 mm OD Cone Tipped Dual Gauge Needle PK6	E	
88000	7105KH 5.0 µL SYR (24/2.75 "/3)	000802	5.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle	E	
88004	7105KHCH 5.0 µL SYR (24/2.75 "/3)	000852	5.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	E	
88007	7105KHWG 5.0 µL SYR (24/2.75 "/3)	000852	5.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Cone Tipped Needle & Repeating Adaptor	E	
88011	7105KH 5.0 µL SYR (24/2.75 "/2)	000803	5.0 µL NanoVolume Syringe with 7 cm 0.63 mm OD Bevel Tipped Needle	E	
88035	75ASN/PE SYR 5 µL 0.63 mm OD	001954	5 µL Fixed Needle PerkinElmer Syringe with 7 cm 0.63 mm OD Cone Tipped Needle	E	
88040	75ASN/PE SYR 5 µL 0.47 mm OD	001953	5 µL Fixed Needle PerkinElmer Syringe with 7 cm 0.47 mm OD Cone Tipped Needle	E	

Every effort has been made to ensure the SGE branded Equivalents and similar close Options suggested here are suitable alternative products, however the appropriateness for individual applications must remain the responsibility of the user.

For more information visit www.sge.com or contact techsupport@sge.com.

