

Syringes | HPLC Autosampler

Diamond Syringe - Syringe Brilliance

Known for 50 years in chromatography syringes, the SGE range of HPLC autosampler syringes are designed and tested to meet critical autosampler specifications.

Product Specifications

Available for CTC, Spark Holland, Hitachi, Thermo Scientific, PerkinElmer, Waters WISP, HPLC autosamplers.

Accuracy and Reproducibility	± 1 % (displaced volume)
Temperature Range	5-70 °C Fixed Needle, 5-120 °C Removable Needle
International Standards Traceability	

For more information about this product visit www.sge.com or contact techsupport@sge.com.

