

# Separations | HPLC Column ProteCol® Cyano

## Turn Imagination into Results

High quality phases and a variety of hardware options provide a perfect match for all your separation needs.

### Features and Benefits

Provides chromatographic retention both in normal and reversed phase separation due to its moderate polarity.

### Recommended Applications

Separation of flavonoids, extraction of polar compounds from non-polar samples as well as the analysis of samples containing analytes with a wide range of hydrophobicity.

### Product Specifications

2.1, 3, and 4.6 mm ID columns available in 50, 100, 150, 250 mm lengths.

Phase	Pore Size (Å)	Particle Size (µm)	Surface Area	% Carbon Load	USP Code
Cyano	120	5	300	5	L10

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).

