

Separations | GC Column BP20 (WAX)


More than five decades of innovative phase technologies and unique fused silica production capabilities, together provide end-to-end separation solutions for all applications.

Features and Benefits

- Industry standard Wax column.
- Bonded and cross-linked.
- Able to be solvent rinsed.

Recommended Applications

- Alcohol.
- Free acids.
- FAME's.
- Aromatics.
- Solvents.
- Essential oils.



Separations | GC Columns

Product Specifications

Polyethylene Glycol phase. Suitable replacement for: DB-Wax, HP-20M, Supelcowax 10, CB-Wax, Stabilwax, Carbowax, HP-Innowax, Rtx-WAX, PE-WAX, RH-WAX, ZB-WAX, TRWAX.



Film Thickness	Operating Temperature
0.1-1.0 μm	20 $^{\circ}\text{C}$ to 260/280 $^{\circ}\text{C}$
> 1.0 μm	20 $^{\circ}\text{C}$ to 240/260 $^{\circ}\text{C}$

For more information about this product visit www.sge.com or contact techsupport@sge.com.