

Separations | HPLC Column ProteCol® Ultra

Turn Imagination into Results

Complementing the ProteCol® analytical range. These columns use 1.8 µm particles for UHPLC use.

Features and Benefits

Specifically designed for UHPLC use, sub 2 micron phases enable use up to 19000 psi.

Recommended Applications

UHPLC.

Product Specifications

The following phases are available in 2 mm ID columns, with 50, 100, and 150 mm lengths:

- ProteCol Ultra C18.
- ProteCol Ultra C8.
- ProteCol Ultra Amino.
- ProteCol Ultra Cyano.
- ProteCol Ultra HILIC FL.
- ProteCol Ultra HILIC PI.
- ProteCol Ultra PFP.
- ProteCol Ultra Phenyl.
- ProteCol Ultra Phenyl Hexyl.
- ProteCol Ultra Polar.
- ProteCol Ultra SCX.

For more information about this product visit www.sge.com or contact techsupport@sge.com.

