

# Connections | GC Ferrules SilTite® FingerTite

The Smart Alternative. One Nut, One Ferrule, All Systems.

SilTite® FingerTite is a GC ferrule system that has been designed to be installed only with the force of you fingers, delivering an easy, leak and maintenance free seal.

## Features and Benefits

- Simplified installation process saves time.
- Leak free - no column creep.
- Cannot overtighten - no broken columns.
- Chemically inert - no contamination.
- Metal ferrule - no temperature limit.

## Recommended Applications

MS, high throughput laboratories.



---

## Product Specifications

Kits available for Agilent, PerkinElmer, Shimadzu, Thermo Scientific, and Varian/Bruker GC-MS systems.

Ferrules available with 0.4, 0.5, and 0.7 mm hole sizes. Blanking ferrules, replacement female nuts, and replacement Agilent base seals also available.

For more information about this product visit [www.sge.com](http://www.sge.com) or contact [techsupport@sge.com](mailto:techsupport@sge.com).